

Certificate in relation to air transportation by passengers with walking disabilities

We hereby certify that the our products

**viamobil V25
e-fix E35 and E36**

including the Lithium ion batteries fulfil the requirements of

IATA chapter 2.3.2.4

and therefore may be carried as checked baggage without restriction.
The lithium ion batteries are tested according UN Manual of Tests and Criteria"
Section 38.3.

Condition precedent to meet the above listed requirement is the transportation in a way that an unintentionally switching on and damage by other luggage is avoided.
The battery must be removed for transportation (protected from damage and from short circuit e.g. by placing the battery in a protective pouch) and **carried** in the **passenger cabin**.

The pilot in command has to be informed about the location of the battery.

We strongly recommend obtaining the approval by corresponding airline company before commencement date respectively departure date.

This certificate does not apply to damaged or defect batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31st December 2019**.

Albstadt, 10th of December 2018

Alber GmbH



Hansjörg Reiner
Managing Director