



Product Service

# CERTIFICATE

No. Z1 13 02 42097 053

**Holder of Certificate:** Invacare Rea AB  
**Org. No. 556098-9211**

Växjövägen 303  
 343 71 Diö  
 SWEDEN

**Production Facility(ies):**

42097

**Certification Mark:**



**Product:** Medical Beds

**Model(s):** Octave

**Parameters:**

Application environment	3, 4
Safe working load	385 kg
Maximum patient weight	350 kg
Dimensions of mattress support platform	119 x 200 cm
Range of height adjustment	40 – 80 cm
Mattress support platform	4-sectioned
Voltage supply	230 VAC / 50 Hz
Max. Current input	3 A
Operating time	10%, max. 2min/18min
Type of applied part	Type B
Protection class	II
Degree of protection	IPX6 / IPX4
Angle A	57°
Angle B	30°
Angle C	115°
Angle D	11°

**Tested according to:** EN 60601-2-52:2010  
 EN 60601-1:2006

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition the certification holder must not transfer the certificate to third parties. See also notes overleaf.

**Test report no.:** 713009210

**Date,** 2013-03-19

(Torsten Zimmer)

Page 1 of 1

