



Product Service

CERTIFICATE

No. B 17 04 11048 025

Holder of Certificate: Alber GmbH

Vor dem Weißen Stein 21
72461 Albstadt
GERMANY

Production Facility(ies): 11048

Certification Mark:



Product: Components for wheelchairs
Push and brake aid for manual wheelchairs

Model(s): viamobil V25

Parameters:

total weight:	13,8 kg
slope:	max. 15,8% at 100 kg user weight max. 11,5% at 160 kg user weight
speed	
- forward:	max. 5,5 km/h
- backward:	max. 2,7 km/h
max. user weight:	160 kg
rated voltage:	36 VDC
power of the motor:	150 W
batteries:	36 V / 6,6 Ah
battery charger:	rated output voltage: 45 VDC rated output current: 1.5 A rated input voltage: 100 - 240 VAC / 50 - 60 Hz input current: 1.35 - 0.7 A

Tested according to: EN 12184:2014

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition the certification holder must not transfer the certificate to third parties. See also notes overleaf.

Test report no.: 713086246

Valid until: 2022-04-03

Date, 2017-04-04

(Jens Meyer)

Page 1 of 1

