



## **Certificate in relation to air transportation of Mobility Aids with Lithium Batteries by passengers with walking disabilities**

We hereby certify that the our products

**viamobil V25  
e-fix E35 and E36**

including the Lithium ion batteries fulfil the requirements of

**IATA chapter 2.3.2.4**

and therefore the system may be carried as checked baggage without restriction.

**Only the battery has to be removed for transportation** (protected from unintentionally switching on, damage and from short circuit by other luggage e.g. by placing the battery in a protective pouch) and **carried** in the **passenger cabin**.

The lithium ion batteries are tested according UN Manual of Tests and Criteria" Section 38.3.

The rated energy of the batterie is **below** the IATA limit of **300 Wh**.

The pilot in command has to be informed about the location of the battery.

**We strongly recommend to contact the airline sufficiently in advance to request approval for transport.**

This certificate does not apply to damaged or defect batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31<sup>st</sup> December 2021**.

Albstadt, 20<sup>nd</sup> of November 2020

Alber GmbH

A handwritten signature in black ink, appearing to be "H. Reiner", written over a horizontal line.

Hansjörg Reiner  
Managing Director