

Certificate in relation to air transportation of Mobility Aids with Lithium Batteries by passengers with walking disabilities

We hereby certify that our products

**viamobil V25
e-fix E35 and E36**

including the lithium-ion batteries fulfil the requirements of

IATA chapter 2.3.2.4

and therefore, the system may be carried as checked baggage without restriction.

Only the battery has to be removed for transportation (protected from unintentionally switching on, damage and from short circuit by other luggage e.g. by placing the battery in a protective pouch) and **carried** in the **passenger cabin**.

The lithium-ion batteries are tested according to UN Manual of Tests and Criteria" Section 38.3.

The rated energy of the batterie is **below** the IATA limit of **300 Wh**.

The pilot in command has to be informed about the location of the battery.

We strongly recommend to contact the airline sufficiently in advance to request approval for transport.

This certificate does not apply to damaged or defect batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31st December 2023**.

Albstadt, 02nd of December 2022

Alber GmbH



Joost Beltman
Managing Director