


# Vicair celleoversigt - gram Smartcells (sølvgrå)

Dato		Fornavn	
Pudetype		Efternavn	
Pudestr.		CPR-nr	


## Vicair Vector O2 9 kamre, 6 cm

Str.	HMI-nr.	1+2	3+4	5/6	7/8	9
3540	105530	48	34	16	39	38
3545	113637	42	29	13	34	33
4040	105537	56	39	18	46	44
4045	105536	62	44	20	51	50
4050	105535	69	48	22	56	55
4343	105534	64	46	21	52	51
4345	113102	52	36	16	42	41
4540	105533	62	44	20	51	50
4545	105532	70	50	22	57	55
4550	105528	78	55	25	64	61
5045	105523	78	55	25	64	61
5050	105531	87	61	27	70	68



## Vicair Vector O2 9 kamre, 10 cm


Str.	HMI-nr.	1+2	3+4	5/6	7/8	9
3540	105519	64	47	26	56	46
3545	113634	55	41	23	48	39
4040	105520	73	54	30	64	52
4045	105521	82	60	34	72	59
4050	105529	91	66	38	80	65
4343	105522	85	63	35	73	60
4345	113093	68	50	28	59	48
4540	105518	82	60	34	72	59
4545	105524	93	68	39	81	65
4550	105525	103	76	43	90	73
4555	113096	87	64	36	76	61
5040	128029	91	67	38	79	64
5045	105525	103	76	43	90	73
5050	105527	115	84	49	99	81
5055	113099	97	81	40	84	68
5543	128032	108	79	45	94	76
5545	128035	113	83	47	98	80
5550	128028	126	92	52	109	89



Yes, you can.®


### Vicair Vector O2 8 kamre, 6 cm

Str.	HMI-nr.	1+2	3+4	5	6+7	8
3540	63443	70	51	34	13	29
4040	63444	80	58	38	15	33
4045	63445	90	65	43	17	38
4050	63452	100	73	48	18	42
4343	63446	93	67	44	17	39
4540	63447	90	65	43	17	38
4545	63448	101	74	49	19	42
4550	63449	113	82	54	21	47
5045	63450	113	82	54	21	47
5050	63451	125	91	60	23	52




### Vicair Vector O2 8 kamre, 10 cm

Str.	HMI-nr.	1+2	3+4	5	6+7	8
3540	61784	89	70	39	21	35
4040	61785	102	80	45	24	40
4045	61786	114	90	51	27	45
4050	61787	127	100	56	30	49
4343	61788	117	92	52	27	46
4540	61789	114	90	51	27	45
4545	61790	129	101	57	30	50
4550	61791	143	112	63	33	56
5045	61792	143	112	63	33	56
5050	61793	159	125	70	37	62



### Vicair Active O2, 10 cm


Str.	HMI-nr.	1 + 4	2 + 3
3540	108714	54	41
4040	108715	62	46
4045	108716	70	52
4050	108717	78	58
4540	108718	70	52
4545	108719	79	59
4550	108720	87	65
5045	108721	87	65
5050	108722	97	73



Yes, you can.®


## Vicair Adjuster O2, 6 cm

Str.	HMI-nr.	For	Midt	Bag
3540	106741	29	77	77
4040	106733	33	89	89
4045	106743	37	100	100
4050	106744	41	111	111
4343	106745	38	102	102
4540	106746	37	100	100
4545	106747	41	112	112
4550	106748	46	125	125
5045	106750	46	125	125
5050	106749	51	138	138




## Vicair Adjuster O2, 10 cm

Str.	HMI-nr.	For	Midt	Bag
3540	106739	42	114	114
4040	106732	48	131	131
4045	106737	54	147	147
4050	106738	60	163	163
4343	106736	56	151	151
4540	106734	54	147	147
4545	106735	61	165	165
4550	106740	68	184	184
5045	106742	68	184	184
5050	106751	75	204	204



## Vicair Multifunctional O2 rygpute


Str.	HMI-nr.	1+9	2+3+7+8	5	4 +6	10
3754	124838	46	52	45	43	39
4054	124839	50	56	49	46	39
4554	124840	56	63	55	52	39
5054	124841	62	70	61	58	39



Velcroåbning på bagsiden til fyldning i

## Vicair Multifunctional O2, 8 cm

Str.	HMI-nr.	For	Midt	Bag
3752	35856	61	115	129
4052	35857	66	124	141
4552	35858	74	139	157
5052	35859	82	155	174



Yes, you can.®