

INVACARE®
TOP END®
SPORT AND RECREATION PRODUCTS



INVACARE®

TOP END®

SPORT AND RECREATION PRODUCTS

In 1986, Invacare Top End was born when founder and designer Chris Peterson began building racing chairs in a small machine shop in Florida. Top End was acquired by Invacare in 1995. The company is now recognized worldwide as the leading manufacturer of high performance mobility products. For more than 28 years, Invacare Top End's mission has remained the same – to be the innovator in sports, recreation and everyday equipment by combining cutting edge materials, engineering advancement and the will to push the boundaries. Invacare Top End's athletes and customers help to achieve this mission by demanding stronger, lighter, faster and more efficient products.

Invacare Top End innovation is the result of a desire to outperform the limits that seemed insurmountable. Innovation starts and ends with the frame, the core of the wheelchair. While significant emphasis is often placed on material choices, the true foundation of the frame is its geometry. Invacare Top End has spent over 28 years perfecting the mathematics that translate into performance and these wheelchairs have a ride and handling that is absolutely legendary.

INVACARE®
TOP END®
SPORT AND RECREATION PRODUCTS

Invacare® Top End® Eliminator™ OSR Racing Wheelchair	4
Invacare® Top End® Crossfire™ All Terrain Wheelchair	6
Invacare® Top End® Pro™ -2 All Sport Wheelchair	8
Invacare® Top End® Pro™ Basketball Wheelchair	10
Invacare® Top End® Schulte 7000 Series Basketball Wheelchair	12
Invacare® Top End® Pro™ Tennis Wheelchair	14
Invacare® Top End® T-5 7000 Series Tennis Wheelchair	16
Invacare® Top End® XLT Handcycle	18
Invacare® Top End® XLT Jr. Handcycle	20
Invacare® Top End® Excelerator™ Series Handcycle	22
Invacare® Top End® Force™ CC Handcycle	26
Invacare® Top End® Force™ -3 Handcycle	28
Invacare® Top End® Force™ RX Adjustable Back Handcycle	30
Invacare® Top End® Force™ K Handcycle	34



INVACARE® TOP END®
ELIMINATOR™ OSR
RACING WHEELCHAIR



The **Invacare® Top End® Eliminator™ OSR Racing Wheelchair** is designed for speed. The super-stiff roll-formed horizontal maintube, combined with a relaxed headtube angle for exceptional stability and aerodynamics, make it the choice of champion racers like Ernst van Dyk, Krige Schabort, Kurt Fearnley, Chantal Petitclerc and Edith Hunkeler. So, if you have your sights set on winning, the Eliminator OSR racing wheelchair will help you get to the finish line first.

Features

6061 T6 oval tube custom aluminum frame
Integral wrap-around fenders
4 cage designs: V, Open-V, U or I-cage
Kneeling, amputee or traditional upholstery
700C or 26" high performance rear wheels
Caliper hand brake
Aluminum side panels
Straight axle tube with rear wheel alignment
Aluminum fork
Precision track compensator system

20" front wheel with tubular tire
Rubber-coated 3/8" diameter aluminum rims
Aero bars
Powder coat paint
250 lb. weight capacity



1. I-cage
2. Open-V cage shown with cut-out axle option
3. V cage
4. U cage shown with knee plate
5. Oversized rake for maximum speed and cornering
6. Precision track compensator system

Options

Carbon fiber rear wheels
Corima™ carbon fiber front wheel
Alignment gauge
Wheel bag
Solid seat/kneel plate
Fixed footrest

Cushion
Cordless speedometer
Helmet
18" front wheel



INVACARE® TOP END®
CROSSFIRE™
ALL TERRAIN
WHEELCHAIR



The **Invacare® Top End® Crossfire™ All Terrain Wheelchair** is the monster truck of wheelchairs and is designed to leave some serious tracks. Do you like to hike, fish, camp, hunt, kayak, garden or go to the beach? You can do it all with the Crossfire All Terrain Wheelchair. Built of 6061 T6 Aluminum, it is equipped with enormous 25" x 2.75" knobby wheels and huge 8" x 2" front casters for extreme traction. The Crossfire All Terrain Wheelchair is engineered with mega ground clearance, plus forward and side-to-side stability so that you can conquer the toughest terrain.

Features

6061 T6 Aluminum frame
Adjustable center-of-gravity
Camber 0-9°
Lock down folding back
Adjustable back height
Adjustable tension seat upholstery
Adjustable tension back upholstery

24" or 25" high performance Sun® wheels
Hard anodized handrims
24" or 25" x 2.75" knobby tires
8" x 2" front casters
Adjustable height and angle footplate
Powder coat paint
250 lb. weight capacity

Model #: CRFAT



1. 24" or 25" x 2.75" knobby wheels and tires
2. 8" x 2" front casters
3. Engineered for front and side-to-side stability
4. Center-of-gravity adjustment
5. Anti-tip option
6. Optional push handles, wheel locks and anti-tip bars

Options

Spinergy® Spox wheels
Quad release axles
Titanium axles
Coated pushrims
Titanium handrims
Natural Fit® handrims
Fixed back
Push handles
Carbon fiber sideguards

Removeable sideguards
Folding sideguards
Wheel locks
Positioning straps
Alignment gauge
Footplate risers
2" cushion
Anti-tip bars



INVACARE® TOP END®
PRO™ -2
ALL SPORT WHEELCHAIR



The **Invacare® Top End® Pro™-2 All Sport Wheelchair** is all about sports. If you want to get into the game, this is the chair for you. Just name it: basketball, tennis, softball, table tennis, badminton, floor hockey or any other sport that is played on a court and the Top End Pro-2 Wheelchair can get you started. The Pro-2 All Sport Wheelchair is offered in three seat widths with lots of adjustments which is perfect for clubs and aspiring athletes. Get in the game.

Features

Aluminum frame: choice of 3 seat widths

Adjustable seat depth 13"-17"

Screw on seat upholstery

Adjustable front and rear seat heights

Adjustable clamp-on tubular recessed footrest

Adjustable center of gravity

Removable single swivel anti-tip bar

Camber: 15°

25" High performance wheels with high pressure clincher tires

Quick release axles

Non-folding adjustable back: angle and height

Adjustable tension back upholstery

Aluminum side guards

Cushion

Positioning strap

Rear frame protectors

250 lb. weight capacity

Model #: P2AS



1. Pro-2 All Sport Wheelchair welded bumper/recessed footrest

2. Adjustable seat depth: 13" - 17"

3. Removable single swivel anti-tip bar

4. Adjustable seat height: front and rear

Options

Basketball carrying bag

Click straps


Metal O-Ring positioning belt

Extra seat cushion

Note: Some assembly is required.



INVACARE® TOP END®
PRO™
BASKETBALL WHEELCHAIR



The **Invacare® Top End® Pro™ Basketball Wheelchair** delivers custom-built performance with adjustments that take only minutes. Perfect for competitors still dialing in their game.

Features

Custom aluminum frame

Offensive 4 bend wing

Adjustable front and rear seat heights

Anti-tip/fifth wheel: Adjustable height

Fixed back

Adjustable back height

Adjustable tension back upholstery

Adjustable Center-of-Gravity

Rear seat height adjustment range is 6" and the frame is available in a tall or short variety to accommodate every aspiring athlete

Fixed camber: 15°, 18° or 20°

High performance Sun® wheels

Anodized pushrims

High pressure clincher tires

Smart fork caster system

72mm casters


Adjustable platform footrest

Powder coat paint

250 lb. weight capacity

Quick release lever allows the player to change his or her rear seat height in just minutes

Model #: SPBB



1. Front seat height and Center-of-Gravity adjustments
2. Rear seat height adjustment
3. Quick release lever
4. Single swivel fixed fifth wheel

Options

Clamp-on fixed aluminum side guards

Spinergy® SLX wheels

Sun® Equalizer™ wheels

12 tab pushrims

Kenda® tires

O-Ring positioning belt

Close mount/cut pushrim tab

Click straps

Spoke guards

Frame protectors

Alignment gauge

Double swivel anti-tip bar

Cushion: soft or hard

Folding back



INVACARE® TOP END®
SCHULTE 7000 SERIES
BASKETBALL WHEELCHAIR




The **Invacare® Top End® Schulte 7000 Series Basketball Wheelchair** is faster, lighter, stronger and better than ever. In 2011 it hit the market with new material, enhanced headtube design and axle securing system, reducing overall weight by two pounds and improving performance. Each Schulte 7000 Series Wheelchair is custom built to each individual athlete's specifications and is overseen by Paul Schulte who is the design engineer for Invacare Top End and the Captain of the USA National Men's Wheelchair Basketball Team.

Features

Custom-built fixed frame in 7005 aluminum
Adjustable height fixed swivel anti-tip 58mm fifth wheel
Fixed back with fixed or adjustable height
Adjustable tension seat and back upholstery
Fixed axle position with fixed wheelbase
Fixed camber: 15°, 18° or 20°
High performance Sun® wheels: 24", 25", 26" or 700C
Four bed offensive wing
Anodized pushrims

High pressure clincher tires
Precision fork caster system
72mm casters
Adjustable height tubular footrest
Powder coat paint
Rear frame protector
Axle securing system
250 lb. weight capacity

Model #: PS7



1. Double swivel anti-tip/ fifth wheel option
2. Precision fork caster system
3. Optional welded fixed tube sideguard
4. Adjustable center-of-gravity option
5. Adjustable height fixed swivel anti-tip
6. 7005 Aluminum heat treating process

Options

Adjustable Center-of-Gravity
Spinergy® SLX wheels
Sun® Equalizer™ wheels
12 tab pushrims
Tapered seat width
Welded seat rigidizer bars with solid seat
Double swivel anti-tip/fifth wheel
Click straps
Double amputee front-end

Recessed platform footrest
Adjustable angle footrest
Clamp-on fixed aluminum side guards
Fixed side guards with welded tube
Fixed aluminum or carbon fiber side guards
Titanium axles
Titanium pushrims
Front frame protectors
Alignment gauge



INVACARE® TOP END® PRO™ TENNIS WHEELCHAIR




For new players, adjustability is key. The **Invacare® Top End® Pro™ Tennis Wheelchair** has just that, featuring adjustable front and rear seat heights, center-of-gravity, footrest angle/height, backrest angle/height and anti-tip. Plus, the Invacare Top End Pro Tennis Wheelchair is available in two frame height choices: short or tall. The short frame is ideal for children and the tall frame is for adults, which makes this chair a great choice for programs that have multiple users. Adjusts in minutes, super easy to use, lightweight and durable. This chair is a winner.

Features

Adjustable front and rear seat heights
Anti-tip/5th wheel: Adjustable height
Fixed back
Adjustable back height
Adjustable tension back upholstery
Adjustable Center-of-Gravity
Fixed camber 15°, 18° or 20°
High performance Sun® wheels:
24", 25", 26"

Anodized pushrims
High pressure clincher tires
Smart fork caster system
72mm casters
Adjustable platform footrest
Powder coat paint
250 lb. weight capacity

Model #: SPTN



1. Front seat height and Center of Gravity adjustments
2. Rear seat height adjustment
3. Quick release lever
4. Single swivel fixed fifth wheel

Options

Clamp-on fixed aluminum side guards
Spinergy® SLX wheels
Sun® Equalizer™ wheels
12 tab pushrims
Kenda® tires
Click straps
O-ring positioning belt

Spoke guards
Frame protectors
Alignment gauge
Telescoping anti-tip bar
Cushion
Folding Back
Fold-in side guards



INVACARE® TOP END®
T-5 7000 SERIES
TENNIS WHEELCHAIR




The **Invacare® Top End® T-5 7000 Series Tennis Wheelchair** is faster, lighter, stronger and better than ever. In 2011 it hit the market with new material, enhanced headtube design and axle securing system, reducing overall weight by two pounds and improving performance. Each Top End T-5 7000 Series Tennis Chair is custom built to each individual athlete's specifications and is overseen by Paul Schulte who is the design engineer for Top End and the Captain of the USA National Men's Wheelchair Basketball Team.

Features

Custom-built fixed frame in 7005 aluminum
Adjustable height fixed swivel anti-tip
58mm fifth wheel
Fixed back with fixed or adjustable height
Adjustable tension seat and back upholstery
Fixed axle position with fixed wheelbase
Fixed camber: 15°, 18° or 20°
High performance Sun® wheels:
24", 25", 26" or 700C
Four bed offensive wing

Anodized pushrims
High pressure clincher tires
Precision fork caster system
72mm casters
Adjustable height tubular footrest
Powder coat paint
Rear frame protector
Axle securing system
250 lb. weight capacity

Model #: T57



1. Custom front frame angle and tubular footrest
2. Adjustable height fixed swivel anti-tip
3. 7005 aluminum heat treating process
4. Precision fork caster system
5. Optional welded fixed tube sideguards
6. Adjustable center-of-gravity option

Options

Spinergy® SLX wheels
Sun Fusion wheels
Adjustable Center-of-Gravity
Welded seat rigidizer bars with solid seat
Telescoping removable anti-tip bar
Click straps
Bolt-on fixed side guards
Fixed side guards with welded tube
Fixed aluminum or carbon fiber side guards
Titanium axles
Titanium pushrims

Coated pushrims
Footrest covers
Frame protectors
Custom strapping/toe clips
Alignment gauge
Adjustable angle footplate/toe clips
Recessed platform footplate with toe stops/footrest strap
Kenda® tires
Close mount/Cut pushrim tab
O-ring positioning belt



INVACARE® TOP END®

XLT HANDCYCLE



If recreation is your goal, the **Invacare® Top End® XLT Handcycle** is great fun and features low maintenance 7 speed hubs with reverse hands-on braking. Just get out there and go! All Invacare Top End recreation handcycles are equipped with a safety flag, full composite chain guard, cushion and seat/chest restraint, as standard, for extra value, comfort and safety.

The **Invacare® Top End® XLT Handcycle** is designed with the recreational handcycling enthusiast in mind. Equipped with seven speeds, the XLT handcycle will easily handle most terrain, whether it's riding around town or taking a spin in the country. The lightweight, oval tube aluminum XLT frame with chrome-moly steel fork features simple geometry for maximum stiffness and easy cornering. The seat and back are totally adjustable for a comfortable fit. The narrow rear-end enhances aerodynamics and the 15° of camber increases stability. If you want a great way to exercise, cross-train or just have fun, choose the XLT handcycle.

Model #: XLT



1. Twist Shifter
2. Quad Twist shifting



Features

Choice of 15" wide or 17" wide seat widths
6061 Aluminum frame/chrome-moly fork
Adjustable seat, back angle and crank height
Full chain guard, chain tension idler and safety flag
Adjustable footrest system with safety straps
Seat and back cushions with chest restraint
Shimano® 7 speed hub with reversing drum brake
175 mm crankset with one chain ring
Quad Twist shifter option helps those with hand impairments

Frictionless ovalized ergonomic hand pedals
Parking brake
26" High Performance wheels
High pressure clincher tires
Quick release axles
15° camber
Road crown compensator
Reflectors
250 lb. weight capacity
Many options and aftermarket accessories are available



INVACARE® TOP END®
XLT JR.
HANDCYCLE



The **Invacare® Top End® XLT Jr. Handcycle** is built just for kids or small adults with 20" wheels and 7 speeds. The XLT Jr. Handcycle's low center of gravity make it stable and a great choice for kids who push the speed limit.

The **Invacare® Top End® XLT Jr. Handcycle** is a recumbent handcycle built just for kids. Easy to operate, this handcycle has seven speeds and an aluminum frame with a chrome-moly steel fork. The enhanced aerodynamics and the 15° of camber increase stability. A tall, narrow back offers maximum range of motion and support, with a chest strap for safety. This extremely stable and maneuverable three-wheel handcycle is great fun for anyone who's between 38" to 60" tall. Standard safety features include: reflectors, full chain guard, footrest safety straps, safety flag, cushion and seat restraint. The XLT Jr. Handcycle is totally adjustable and ready to ride!

Model #: XLTJR



1. Parking brake
2. Optional horizontal handpedals
3. Ergonomic handpedals



Features

- | | |
|---|--|
| 6061 Aluminum frame/chrome-moly fork | Parking brake |
| 14" wide adjustable seat, back angle and crank height | 20" cruiser wheels and tires |
| Full chain guard | Quick release axles |
| Chain tension idler and safety flag | 15° camber |
| Adjustable footrest system with safety straps | Road crown compensator |
| Seat and back cushions with chest restraint | Reflectors |
| Shimano® 7 speed hub with reversing drum brake | 250 lb. weight capacity |
| 155 mm crankset with one chain ring | Many options accessories are available |
| Vertical ergonomic handpedals | |



INVACARE® TOP END®
EXCELERATOR™ SERIES
HANDCYCLE



If you want a great way to exercise, cross-train or just have fun, one of the **Invacare® Top End® Excelerator™ Series Handcycles** is just what you need. This upright and extremely maneuverable three-wheel handcycle is capable of speeds up to 15 mph. The smooth seven-speed Shimano® hub with its reversing drum brake delivers hands-on braking and allows the rider to tackle hills or level terrain with equal ease. The oval tube carbon steel frame construction with adjustable sliding seat makes transfers in and out easier too.

The Excelerator Handcycle is built to your specifications, plus you can order the popular mountain drive option which adds 7 extra low gears for climbing hills or the V crankset option with its wider cranks to allow more clearance for those who are tall or require additional space. No excuses...Get Cranking!

Invacare® Top End® Excelerator™ Features

Carbon steel frame and fork
Adjustable seat: forwards and backwards
Adjustable back height
Fold over back upholstery
Adjustable crank height
Chain tension idler
Full chain guard, reflectors and safety flag
Adjustable height/angle footrest system
Shimano® 7 speed hub w/reversing drum brake

Parking brake
175 mm crankset w/one chain ring
Ovalized ergonomic hand pedals
Performance wheels w/cruiser tires
Quick release axles
3° camber
Cushion and hook and loop seat restraint
Weight capacity 350 lb.

Model #: XCL



1. Twist Shifter for 7 speed hub
2. Quad Twist shifter option
3. Parking brake
4. Mountain drive transmission option
5. Frictionless ovalized handpedals



Options

Mountain drive transmission adds 7 low gears
Quad twist shifter helps those with hand impairments
Tri-pin quad handpedals
Crank width adapter set
Quad gloves
Water bottle and cage
Tow bar for wheelchair
Cordless computer

Front safety light
Rear flashing light
Helmet
Seat restraint: Auto style fastener
Mirror
Crutch holder and strap
Bike rack for trailer hitch
Top End t-shirt

Invacare® Top End® Li'l Excelerator™ Features

Carbon steel frame and fork

Adjustable seat: forwards and backwards

Adjustable back height

Fold over back upholstery

Adjustable crank height

Chain tension idler

Full chain guard, reflectors and safety flag

Adjustable height/angle footrest system

Shimano® 7 speed hub with reversing drum brake

Parking brake

155 mm crankset with one chain ring

Foam covered horizontal handpedals

20" Cruiser wheels and tires

Quick release axles

9° camber

Cushion and hook and loop seat restraint

Weight capacity 250 lb.

Model #: LXC



1. Parking brake

2. Ergonomic handpedals

3. Hook and loop safety seat restraint



Options

Quad twist shifter option

Quad gloves

Mirror

Water bottle and cage

Helmet

Front safety light

Rear flashing light

Crutch holder and strap

Bike rack for trailer hitch

Seat restraint: Auto style fastener

Features

1



1. V Crankset

The Invacare® Top End® V crankset option accommodates consumers who might have “clearance” issues. Not compatible with mountain drive option. Available on the Invacare Top End Custom Exceleator Handcycle.

2



2. Mountain Drive Transmission

Available on the Exceleator handcycle
Adds 7 low gears (7 low + 7 standard = 14 total)
Easy shifting, especially for quadriplegics – just push-in the center button
Lightweight planetary transmission
A must have for hilly terrain

3



3. Full Chain Guard

A full chain guard on all 7 speed handcycles is a great safety feature



INVACARE® TOP END® FORCE™ CC HANDCYCLE



The **Invacare® Top End® Force™ CC Handcycle** is a handcycle you can ride on wide trails that is lightweight but also includes extreme climbing gears, sturdy 559-26" bike wheels equipped with mountain bike tires and a stop-on-a-dime front disk brake. The design is a cross between the popular Invacare® Top End® Force™ RX Handcycle and the Invacare® Top End® Force™ G Handcycle with enhanced ground clearance and an aggressive back angle adjustment that puts you in charge of rough terrain. Now you can force the pace....cross country especially with the BionX® Power Boost option which adds 25%, 50%, 75% or 100% power assist.


Features

Ovalized 7000 series aluminum frame
Front disk brake
Adjustable seat length
Adjustable back angle
Aluminum fork
Adjustable crank housing
Multi-position adjustable footrests
Safety straps
Invacare® Top End® V crankset

Ovalized aluminum ergonomic handpedals
Shimano® 27 speed drive train
Shimano® Rapid Fire™ hands-on shifter brake
11-34 SRAM® cassette
Cable driven upper derailleur shifter
Upper chainrings: 44-34-22
Steering dampener
559-26" high performance wheels

Knobby or Cruiser tires
Threaded axles
Parking brake
Seat and back cushions
Chain guard
Safety flag
250 lb. weight capacity
Stealth Black powder coat paint
BionX® Power Boost option adds 25%, 50%, 75% or 100% power assist

Model #: FRCCC



1. Front disk brakes
 2. BionX® Power Boost (available in the US and Canada only)
 3. Optional dual rear disk brakes
 4. Hands on shifting and braking
 5. Welded footrest clamps
 6. Brake lever for optional rear disk brakes
-

Options

Heart rate monitor
Alignment gauge
Backpack hydration system
Rear disk brakes
Water bottle and cage
Cordless computer
Helmet
Safety lights, mirror
Additional safety flag
Click straps

Bike rack
Indoor training roller
Extra cable kit
Leg guard attachment
Tool and tire repair kit
Wheel bag
Handcycle travel bag
Top End Bike Jersey
Crank width adapter set



INVACARE® TOP END®
FORCE™ -3
HANDCYCLE



The **Invacare® Top End® Force™ -3 Handcycle** is a great handcycle at a great price. This bike is ready to ride with 27 speeds, performance wheels and Shimano® components. Perfect whether you are just looking to get started in handcycling or ready to force the pace at your next competition.

Features

15" or 17" Width

Sealth Black Power Coat Paint

60" Wheelbase

Adjustable back angle, seat and footrests

27 speed drive train

Shimano® Deore Rapid Fire™ shifter/brake

Shimano® Octalink 170 mm cranks

Shimano® Tiagra brakes

Shimano® Sora manual upper derailleur

Adjustable crank height: diagonal and verticle

Alex DM18 silver spokes 26-559 wheels
9° camber

High pressure cruiser clincher tires

Seat and back cushions

Chest and footrest safety straps

Drafting bumper

Threaded axles

Frictionless ovalized ergonomic hand pedals

Safety flag

Chain guard


Leg guard

Reflectors

Parking brake

250 lb. weight capacity

Model #: FRC3



1. Wrap around frame with adjustable seat upholstery
2. Cable driven shifting for upper chainrings
3. Front disk brake and welded footrest clamps
4. Backrest adjusts forwards and backwards
5. Adjustable back angle allows 35-55° recline
6. Shimano components for hands-on shifting and braking

Options

Heart rate monitor

Alignment gauge

Backpack hydration system

Water bottle and cage

Cordless computer

Helmet

Safety lights, mirror

Crank width adapter set

Additional safety flag

Click straps

Bike rack

Indoor training roller

Extra cable kit

Tool and tire repair kit

Wheel bag

Top End Bike Jersey



INVACARE® TOP END® FORCE™ RX HANDCYCLE



The **Invacare® Top End® Force™ RX Handcycle** with adjustable carbon fiber back is the latest competitive handcycle in the Invacare® Top End® Force™ Series. The Force RX Handcycle combines the hypersonic speed of the Invacare® Top End® Force™ X Handcycle design with the adjustability of the Invacare® Top End® Force™ R Handcycle in a lightweight aerodynamic package. The 7005 aluminum heat treated frame is fabricated with slick one piece side frames, plus there are multiple enhancements to allow a more reclined position without interferences. The custom reclined frame with adjustable feather weight carbon fiber back is ultra-light and super stiff. The Force RX Handcycle is a completely customized racing machine which includes a SRAM®X9 30 speed drive train, or you can upgrade to top-of-the-line SRAM XX or Rival Shifters. Are you ready to force the pace? Get the new Invacare Top End Force RX Handcycle!


Features

7005 aluminum aerodynamic frame
Adjustable Carbon fiber back
7005 aluminum aerodynamic fork
Adjustable height crank housing
Breathable mesh upholstery
Seat and back cushions
Multi-position adjustable footrest
Tapered aluminum crank arms
Frictionless pedal bearings
Aluminum ergonomic hand pedals

SRAM® X9 30 speed drive train
Chainrings: 48-39-24
Cassettes: SRAM 11-34
Cane Creek SCR-3 Caliper Brake
Fork mounted locking second brake
Steering dampener
26" high performance wheels
Reflectors
High pressure clincher tires
0° or 3° of camber

Threaded axles
Drafting bumper
Parking brake
Positioning straps
Safety flag
Carbon fiber chain guard
9 paint choices
Reflective stickers
250 lb. weight capacity

Model # FRCRXA



1. Derailleur post non-locking second brake is a no charge option
2. Sleek new drop outs are slim, lightweight and aerodynamic
3. Carbon Fiber chain guard looks great and keeps you safe

Upgrades and Options

Carbon fiber wheels
20" rear wheel option
Component upgrades: SRAM XX or Rival Shifters
Fork mounted or handle mounted
SRAM Cassette: 11-28
Chainrings: 53-39-28 or Compact 110
Disk brakes
Quad cuff power plates
Double Tap elbow shifting and braking system

Continental® Grand Prix high pressure clincher tires
Continental Sprinter tubular tires
Bike Jersey
SRAM® Training system
Travel bags
Safety options and hydration systems

See order form for complete details

Features



1. The totally adjustable featherweight carbon fiber back is as lightweight as the fixed welded back, allowing athletes the ability to find their most aerodynamic and powerful position for the course, whether it is a road race, time trial or triathlon.

The bolt-on adjustable headrest cradles the cervical spine and is adjustable in height and angle for a perfect fit. The carbon fiber back is covered with a breathable mesh material that is both lightweight and durable.



2. Choose from 3 fork styles:

The fixed length/vertical height adjustment fork is standard. This fork features the best of both worlds – a super stiff fork and height adjustments so you can get the clearance you need and keep the crank height as aerodynamic as possible.

The fixed height fork is optional and the stiffest, but it is not adjustable in any way. Measurements are required, so it is best to have this dialed in before you go for this style (shown).

The diagonal and vertical height adjustment fork is optional and will give you the most adjustability, especially with the welded back, but it is not as stiff as the other two styles.



3. The SRAM® X9 30 speed drive train, trigger shifter and brake mounted on the right handpedal is standard. This allows for hands-on shifting and braking right at your fingertips.

The ovalized aluminum ergonomic hand pedals are available in 3 sizes: small, medium or large perfect fit. They are offered in 15° or 25° angles.

The power plate option (shown) allows full contact through the entire pedal stroke, especially the push phase.

The standard front wheel cassette is SRAM 11-34, perfect for hilly terrain. You can get a SRAM 11-28 which delivers the closer gears preferred by elite athletes, but it is less of a fit for climbing.



4. Component upgrades are also available:

SRAM double tap 10 speed Rival fork mounted shifter and brake (shown). To use this system you need to have your arms in a position "ready" to brake as the brake is not located on the handpedal, but it does eliminate the shifter and brake cables. With this system you can have a second brake mounted on the left handle, or a fork mounted one-handed locking brake or derailleur post non-locking brake (shown in picture 1 on page 32).

SRAM XX 10 speed right handle mounted shifter and Pro Force brake lever deliver hands-on shifting and braking just the same as the SRAM X9.



5. The redesigned curved monocoque front triangle is super strong and looks terrific. This new design allows the athlete to get closer to the front wheel which translates to a lighter and stiffer fork because it is shorter. Plus, it allows for more body clearance so the athlete can get lower and obtain a more streamlined position.



6. Wheel choices: The 26" High Performance 36 spoked wheels with threaded axles are standard and equipped with high performance clincher tires and tubes. (Not shown)

A set of carbon fiber wheels will definitely upgrade your ride because they are stiffer, lighter and more aerodynamic.

The Top End Rim Upgrade (shown) has a carbon composite rim and stainless spokes: 24 front and 20 rear. These wheels require tubular tires.

The Top End/Spinergy Rim Upgrade has a carbon/graphite composite rim and requires tubular tires.

The Time Trial Package includes the 26" Top End Carbon Rim front wheel and a pair of 26" Zipp Discs on the rear. These are for use in the time trial competition only and are not legal for road racing.



7. Options for cyclists with limited hand function are the Double Tap Frame Mounted Shifters which are positioned so that the right shifter operates lower derailleur and the left operates upper chainrings by using your elbows. Add the Elbow Activated Frame Mounted Brake and you are good to stop and go! A disk brake is recommended for added safety and "stoppability."



8. The quad cuff power plate handpedals feature an adjustable quad cuff to accommodate various hand sizes and positions. There is a detachable strap that fastens to the quad cuff if needed. The lower power plate will help you power through your stroke on the downstroke and the upper power plate will do the same on the upstroke. Most athletes customize their handpedals once they start riding by adding bicycle handlebar tape to optimize their pedal power.



INVACARE® TOP END® FORCE™ K HANDCYCLE



The **Invacare® Top End® Force™ K Handcycle** is a kneeling position handcycle available in two custom, made-to-measure aerodynamic designs. The Invacare Top End Force K Kneeler is used by the most able athletes, while the Invacare Top End Force K Double Amputee Kneeler is custom designed for the double above-the-knee amputee. What makes the Invacare Top End Force K Handcycle different from the rest of the Invacare Top End Force Series handcycles? The Invacare Top End Force K Handcycle is a kneeling-only, non-adjustable position handcycle that requires dialed-in measurements and engineering drawings—all of which are easy to find on the website: www.topendwheelchair.com. Can you force the pace? "Yes, you can.®" Get the Invacare Top End Force K Handcycle.

Features

7005 aluminum custom welded frame
7005 aluminum aerodynamic fork
Breathable mesh upholstery
Seat and back cushions
Tapered crank arms
Frictionless pedal bearings
Aluminum ergonomic hand pedals
SRAM® 30 speed drive train
Chainrings: 53-39 or 48-39-24
Cassettes: SRAM® 11-34

Cane Creek SCR-3 Caliper Brake
Fork mounted locking second brake
Steering dampener
26" high performance wheels
0 or 6 degrees of camber
Threaded axles
Drafting bumper
Parking brake
Safety flag
Chain guard

9 paint choices
Reflective stickers
250 lb. weight capacity

Model #: FRCK



1. Top End Carbon Wheel Upgrade is stiff, light and aerodynamic
2. SRAM components allow shifting and braking at your fingertips

Options

Carbon fiber wheels
Component, brake and cassette upgrades
SRM™ Power Meter measuring tool integrated directly into Invacare Top End cranks

See the Force K order form for all of the great options available!

Frame Styles:

Invacare Top End Force K Kneeler designed for the most able athletes with a custom knee/seat bucket.

Invacare Top End Force K Double Amputee Kneeler designed for the double above-the-knee amputee utilizing a custom cage.

Invacare Corporation

www.invacare.com

Top End

www.topendwheelchair.com

DENMARK

Invacare
Sdr. Ringvej 37
2605 Brøndby
tlf 3690 0000 Fax 3690 0001
denmark@invacare.com
www.invacare.dk

©2014 Invacare Corporation. All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Rotor is a registered trademark of Rotor Components, LLC. Corima is a trademark of Corima, S.A. Zipp is a registered trademark of Composittech, Inc. Spinergy is a registered trademark of Spinergy Inc. Natural Fit is a registered trademark of Three Rivers Out Front. Sun and Equalizer are registered trademarks of Sun Components. Kenda is a registered trademark of Kenda Rubber Industrial Co., Ltd. Shimano and Rapid Fire are registered trademarks of Shimano Inc. SRM is a registered trademark of SRM Service Center Inc. USA. SRAM is a registered trademark of SRAM Corp. Continental is a registered trademark of Continental Tire Corp. BionX is a registered trademark of BionX International Corporation. Specifications are subject to change without notification. Form No. PB-600-EU 131219

