



Product Service

CERTIFICATE

No. B 14 06 11048 015

Holder of Certificate: Alber GmbH

Vor dem Weißen Stein 21
72461 Albstadt
GERMANY

Production Facility(ies): 11048

Certification Mark:



Product: Propulsion systems
Additional drive for manual wheelchairs

Model(s): twion M24

Parameters:

rated voltage:	37 VDC
motor rated power:	2 x 60 W
max. speed of the propulsion support:	6 km/h
max. safe slope:	6°
max. transport weight:	120 kg
degree of protection:	IPX4
batteries:	Li-Ion (37 V, 2,25 Ah)
battery charger:	IN: 100-240 VAC, 50-60 Hz; OUT: 0.83 A - 0.36 A

Tested according to: DIN EN 12184:2009
ISO 7176-14:2008

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition the certification holder must not transfer the certificate to third parties. See also notes overleaf.

Test report no.: 713040387

Valid until: 2016-01-10

Date, 2014-06-16

(Torsten Zimmer)



Page 1 of 1